

Fight the Flu!


ThermoStrip Flu Kit with Symptom Chart

Place thermometer on forehead for 15 seconds. Read yellow moving line.



DISPOSABLE

and



REUSABLE

118F ThermoStrip Flu Kit

Thermometers nestled inside folded chart and instruction sheet.

Size: 4" w x 2.5" h Imprint: 3" w x 1.75" h

Weight: 12 lbs./M

Qty/Item	250	500	1000	2500	5000
118F	\$2.40	2.28	2.17	2.00	1.86

Price includes 1 color imprint on back of instructions.

- Self-application
- Minimal patient caregiver contact



KNC
Marketing
CREATIVE SOLUTIONS

301.843.1896 phone • www.knc.com

Prices include one color imprint/one location. Setup, taxes & freight additional.

Offer expires 1.31.10. Use promo code: Therm10 when ordering.